HOMEDICS[®]

THERA·P[®]

NON-CONTACT

INFRARED BODY THERMOMETER

Model: HTD8816C

2 YEAR GUARANTEE

FOREWORD

The Non-Contact Infrared Body Thermometer operating Instructions intend to provide the necessary information for proper operation of the thermometer model.

Only body mode was reviewed and certified by notified body. General knowledge of Infrared thermometer and an understanding of the features and functions of the thermometer model are prerequisites for proper use. The Non-Contact Infrared Body Thermometer is a medical device, and can be used repeatedly for the product's recommended lifetime of 5 years.

Please read the manual first before using it, if not fully understand the usages, please stop using the thermometer.

Do not operate the thermometer without completely reading and understanding these instructions.

Notice

Purchase or possession of this device does not carry any express or implied license to use with replacement parts which would, alone or in combination with this device, fall within the scope of one of the relating patents.

For further information contact:

FKA Brands Ltd, Somerhill Business Park, Five Oak Green Road,
Tonbridge, Kent TN11 0GP.

Customer Service: service@homedics.co.uk

Share Info Consultant Service LLC Repräsentanbüro
Heerdter Lohweg 83, 40549 Düsseldorf

SAFETY INFORMATION

This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application. The Non-Contact Infrared Body Thermometer is designed to minimize the possibility of hazards from errors in the software program by following sound and light engineering design processes, Risk Analysis and Software Validation.

WARNING

Warnings are identified by the WARNING symbol shown above.

- The Non-Contact Infrared Body Thermometer is to be operated by consumers in the home setting and primary care setting as screening tool. This manual, accessories, Directions for Use, all precautionary information, and specifications should be read before use.
- This product is designed to measure human body temperature on the forehead. Do not use it for any other purpose.
- This product is intended for use in the home setting and primary care setting as screening tool.
- Do not use the thermometer if it malfunctions or has been damaged in any matter.
- When the ambient temperature of the thermometer changes too much, such as moving the Thermometer from one place of lower temperature to another place of higher temperature, allow the thermometer to remain in a room for 30 minutes where the temperature is between 15°C to 40°C (59°F- 104°F).
- Remove primary batteries if equipment is not likely to be used for long time.
- This product is not waterproof, do not immerse in water or other liquid; If cleaning and disinfecting, please follow the "Care and Storage" section requirements.
- Do not touch the sensor of infrared detection with your fingers.
- Avoid using the thermometer to take a body temperature reading while using a cold compress or other measure to reduce fever.
- If measuring human forehead temperature, please select "body" mode; for measuring other objects, liquids, food and other temperature please select "surface" mode.
- This product must be operated in a stable environment. If the ambient environment is unstable, please check whether there is fog on the sensor and remove this before using in accordance with the "Care and Storage" section.

- Do not near strong electrostatic field or strong magnetic fields, thus avoiding the impact on the accuracy of the measurement data.
- Do not mix the old and new batteries to avoid damage to the product.
- It may affect the accuracy of measurements when the forehead is covered by hair, perspiration, cap or scarf.
- The measuring result of this product is only for your reference. If you have any doubt, please measure the temperature in other methods.

⚠ The device should be kept out of the reach of children/pets. When not in use, store the device in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the thermometer body.

⚠ Do not throw batteries into fire.

⚠ Only use recommended batteries. Do not use rechargeable batteries.

⚠ This thermometer does not replace a hospital diagnosis.

⚠ Do not drop, disassemble or modify the device.

⚠ Do not use this device if you think it is damaged or notice anything unusual.

⚠ This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the 'Technical Specifications' section.

⚠ Not servicing/maintenance while the thermometer is in use.

⚠ When using, shall not touch battery and the patient simultaneously.

⚠ Do not use the device if it is damaged/ degraded/loosened in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.

⚠ Based on the current science and technology, other potential allergic reactions are unknown.

⚠ This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS.

OVERVIEW

Intended Use

The HoMedics TheraP Non-Contact Infrared Body Thermometers are designed to be used for intermittent measurement and monitoring of human body temperature by consumers in the home setting and primary care setting as screening tool.

Description of Non-Contact Infrared Body Thermometer

- Device principle and introduction
- The HoMedics TheraP Non-Contact Infrared Body Thermometer is a hand-held, reusable, battery operated device, which can measure human body temperature on forehead, the skin temperature on one's forehead.
- The operation principle is based on Infrared Sensor technology. The IR sensor can output different signal when measuring different object temperature or in different ambient temperature, and the ASIC can turn the signal from IR Sensor to a digital value and display it on the LCD.
- Description on Controls, Indicators, and Symbols


Figure 1: Non-Contact Infrared Thermometer

1. IR sensor
2. Liquid crystal display (LCD)
3. MODE button
4. MEMO button
5. ON/scan button
6. Battery cover
7. Surface mode
8. Body mode
9. Data indicator
10. Indicator of measurement results
11. Low Battery indicator
12. Voice on/off indicator
13. Memory number
14. Fahrenheit
15. Celcius
16. Memory indicator

THERMOMETER APPLICATIONS

| Thermometer Model Number | Thermometer Style | Adult | | Pediatric | |
|--------------------------|--|-------|----------|-----------|----------|
| | | Ear | Forehead | Ear | Forehead |
| HTD8816C | Non-Contact Infrared Body Thermometers | | ✓ | | ✓ |

EQUIPMENT SYMBOLS

| | | | |
|---|---|---|---|
|  | Warning |  | Restriction of Hazardous Substances |
|  | Non-sterile packaging |  | Batch code |
|  | Refer to operating instructions |  | Compliance with WEEE Standard |
|  | Operating Temperature |  | DO NOT THROW AWAY Intended for multiple use |
|  | Operating Humidity |  | Operating atmospheric pressure |
|  | This device complies with Part 15 of FCC (Federal Communications Commission) Rules. |  | Indicates this device is in compliance with MDD 93/42/EEC. 0598 is the Notified Body Number |
|  | Manufacturer |  | Authorized Representative in the European community |
|  | Recyclable |  | Serial number |
| IP22 | IP22: The first number 2: Protected against solid foreign objects of Φ 12.5 mm and greater. The second number: Protected against vertically falling water drops when enclosure tilted up to 15°. | | |

TECHNICAL SPECIFICATIONS

| | |
|--|--|
| Measurement unit | °C/°F |
| Operating mode | Adjusted mode (Body mode) Direct mode (surface mode) |
| Reference Body Site | Axillary |
| Rated output range | Body mode: 34.0 - 43.0°C / 93.2 - 109.4°F Surface mode: 0 - 100.0°C / 32 - 212°F |
| Output range | Body mode: 34.0 - 43.0°C / 93.2 - 109.4°F Surface mode: 0 - 100.0°C / 32 - 212°F |
| Laboratory Accuracy | Body mode: 34.0 - 34.9°C: ± 0.3°C / 93.2 - 94.8°F: ± 0.5°F 35.0 - 42.0°C: ± 0.2°C / 95.0 - 107.6°F: ± 0.4°F 42.1 - 43.0°C: ± 0.3°C / 107.8 - 109.4°F: ± 0.5°F |
| | Surface mode: ± 2°C / ± 3.6°F |
| Display Resolution | 0.1°C/0.1°F |
| Three-colour Backlight (Colour Alarm) | 35.5 - 37.3°C / 95.9 - 99.1°F: Green (Normal Temperature) 37.4 - 38.0°C (Alarm point) / 99.3 - 100.4°F: Yellow (Slight Fever) 38.1 - 43.0°C / 100.6 - 109.4°F: Red (High Fever) Note: 1. Surface mode is always with Green backlight 2. In body mode 34.0-35.4°C is with green backlight |
| Auto Power Off Time | ≤18s |
| Measuring Time | ≤ 2s |
| Measuring Distance | 1 - 5cm (0.4 - 2 in) |
| Memories | 50 |

| Power Supply Requirements | | | |
|---------------------------------|--|-----------------------------|------------------------------------|
| Batteries | 2 x 1.5V (AAA) Alkaline battery (IEC Type LR03) | | |
| Adaptable Range | 2.6V - 3.6V | | |
| Environmental | | | |
| Operating Condition | Operating Temperature: 15 - 40°C (59 - 104°F) | Relative Humidity: ≤ 85% | Atmospheric pressure: 70-106Kpa |
| Transport and Storage Condition | Storage Temperature: -20 - 55°C (-4 - 131°F) | Relative Humidity: ≤ 93% | Atmospheric pressure: 70-106Kpa |

| Dimension and Weighting | |
|--------------------------------|---|
| Weight (without batteries) | 80g |
| Size | L:122mm X W:96mm X H:42mm |
| Compliance | |
| Item | Compliant with |
| Equipment classification | Safety Standards: EN 60601-1: 2006+A1:2013, EN 60601-1-2: 2015 |
| Type of protection | Internally powered equipment (on battery power) |
| Degree of protection | Non-Applied part |
| Front panel and case labelling | EN ISO15223-1:2016 |
| Temperature | EN ISO80601-2-56:2017 |
| Home healthcare environment | EN 60601-1-11:2015 |

Calculated values of the indicators according to ISO 80601-2-56

| | Group A1 | Group A2 | Group B | Group C |
|---|----------|----------|---------|---------|
| Bias | 0.015 | 0.000 | -0.042 | -0.040 |
| Standard deviation | 0.123 | 0.121 | 0.121 | 0.137 |
| Limits of agreement | 0.245 | 0.243 | 0.243 | 0.274 |
| Clinical repeatability (for all groups) | 0.071 | | | |

Note: the above value is calculated from clinical data of HTD8818A.

Safety classification of ME EQUIPMENT

| | |
|---|---------------------------------|
| Protection against electric shock | Internally powered ME equipment |
| Applied part | Non-Applied part |
| Protection against harmful ingress of water or particulate matter | IP22 |
| Mode of operation | Continuous operation |

Note: Not intended to be sterilized. Not for use in an OXYGEN RICH ENVIRONMENT.

OPERATION

Battery installation

Caution: The Non-Contact Infrared Body Thermometer does not operate with dead batteries or external power sources. Install new batteries.

1. Slide the battery door downward, toward the bottom of the Non-Contact Infrared Body Thermometer, and remove the battery access door;
2. Insert two pieces AAA size batteries according to the "+" and "-";
3. Close the battery cover.

How to operate

Before applying the Thermometer, be sure to read and understand all warnings listed of the instructions before use.

- The thermometer is aligned with the middle of the forehead to measure body temperature between the eyebrows above and keep the vertical distance, press the On/scan button, temperature will display immediately, see figure 2.


Figure 2: Measuring position and distance

- When the ambient temperature of the thermometer changes too much, such as moving the Thermometer from one place of lower temperature to another place of higher temperature, allow the thermometer to remain in a room for 30 minutes where the temperature is between 15°C to 40°C.
- The ambient temperature around the test person should be stable, should keep away from the larger flow fan, air-conditioning vents and so on.
- When people moving from one place of lower temperature to another place of higher temperature, should at least remain in the test environment more than 5 minutes, to be consistent with the ambient temperature after the re-measurement.

- Wait at least 1 second for the next measurement. If the continuous measurement of five times, it is recommended to wait at least 30 seconds and then continue measurement.
- You cannot use the thermometer in place where the sun is strong.
- If for some reason the forehead temperature measurement is low, try measurement behind the ears. See figure 3.


Figure 3: Align behind the ears to measurement

GENERAL SET UP AND USE

Start measuring

1. Turn on the thermometer by pressing the On/scan button. The thermometer will perform self-test with all segments displayed for 2 seconds.
2. Align the thermometer with the forehead, keeping the recommended distance, and then press the ON/scan button to start the measurement. Data reading will be displayed.

Note:

1. After measurement you will hear a "beep beep beep" three times which means that the measurement is complete. The temperature is displayed on the LCD. The LCD backlight colour will change red, green or yellow according to the measurement. Green means Normal Temperature (35.5-37.3°C / 95.9-99.1°F) ready for next measurement. Yellow means Slight Fever Warning (37.4-38.0°C / 99.3-100.4°F) Please pay attention to body temperature. Red means High Fever (38.1-43.0°C / 100.6-109.4°F) Please take action to cool down or go to the doctor.
2. To ensure the accuracy of the measurement, wait at least 30 seconds after 5 consecutive measurements.

Mode conversion

- When the device is running, pressing the MODE button to cycle conversion between "body" mode and "surface" mode.
- "Body" mode is used for measuring human body temperature, the "surface" mode is used to measure the surface temperature. (The factory default is "body" mode)

Recalling and Erasing Memory Data

The last temperature taken before the thermometer powers off is stored in memory, up to 50.

1. In the boot or shutdown state, short press the MEMO button to view the history of measured values.
2. An empty memory cell shows "---°C" or "---°F".
3. Temperature readings can be stored in memory. Up to 50 temperature readings can be stored into the memory cells and automatically overwrite historical data.
4. In boot mode, press the MEMO button until the LCD display "CLR", which means that all stored data is cleared completely.

Parameter settings

As below step:

1. Unit Set-F1

Under the boot mode, long press the "MODE" button to enter F1, press the "MODE" or "MEMO" button to switch Celsius and Fahrenheit temperature units, press the ON/Scan button to confirm the unit settings (factory default is Celsius).

2. Fever alert set-F2

Under F1 state, press ON/scan button to enter the F2, press the "MEMO" button to decrease 0.1°C, press the "MODE" button plus 0.1°C, long press to accelerate the speed of temperature regulation, and finally press the ON/scan button to save. (The factory default is 38.1°C)

3. Prompt sound settings-F3

Under F2 state, short press ON/scan button to enter F3, press MODE button or MEMO button to set voice switch, and press the ON/scan button to confirm the settings. (The factory default is the voice Prompt to open).

4. Exit setting mode

In the F3 mode, press the On/Scan button will automatically turn off the screen, exit setting.

- **Prompt sound settings ON/OFF function**

Under the boot mode. Short press on/off button, at the same time press the "MODE" button to take on or take off the Prompt sound function.

- **Restore to factory setting function**

Under the boot mode, short press on/off button, at the same time press the MEMO button until LCD display "rst". Two seconds later, former F1-F3 parameter back to factory setting.

TROUBLESHOOTING

| Message | Situation | Solution |
|---|--|---|
|  Hi | Temperature taken in not within Typical human temperature range. (34.0°C - 43.0°C or 93.2°F - 109.4°F). | Make sure the forehead thermometer is for forehead measurement, not other human body site. |
|  Lo | | |
|  Lo | Measured over the distance: 1-5cm (0.4-2in) | Optimum measurement distance is 1cm |
| | Incorrect test position. | See figure 2 Measuring position and distance. |
| | Subjects forehead hair, Antipyretic stickers, head with sweat, etc. | Subjects sit quietly 5-10 minutes before the test. |
| | Some people's body temperature is lower than the general population. | The main concern fever temperature |
|  Err | Operating temperature exceeds the range of specified temperature. | Move to a room within the operating range wait 30 minutes before taking temperature. |
|  | The screen flicker, automatic turn off. | Replace battery or the product has been damaged, needs repairs. |
|  | Battery capacity is too low. Taking Temperature is not allowed. | Install a new battery |
|  POS | Ambient temperature changes too fast | Wait until the ambient temperature is stable. |
|  | <ol style="list-style-type: none"> Power is off. Improper battery installation. The battery is exhausted. Display remains blank. | <ol style="list-style-type: none"> Press ON button again. Check the battery polarity. Replace with a new battery. Contact the retailer or service center. |

REPLACING THE BATTERY

- Open and release battery cover following indicator on the surface of battery cover. Before changing the battery be sure the thermometer is already turned off.
- Remove the battery and replace with 2 new type AAA batteries, Make sure to align properly as indicated inside the battery cover.
- Slide the battery cover back in until it snaps in place.
Do not dispose of used batteries in household waste. Take them to special local collection sites.

4. If the thermometer appears frozen after changing the batteries, remove batteries, wait 30 seconds and replace batteries again.

WARNING

Do not recharge, disassemble or dispose of in fire.

1. The typical service life of the new and unused batteries is 2000 measurements for the operation time is 18s.
2. Only use the recommended batteries, do not recharge non-rechargeable batteries and do not burn them.
3. Remove the batteries if the thermometer is not to be used for a long period.

CLEANING, CARE AND STORAGE

- The lens is very delicate.
- It is very important to protect the lens from dirt and damage.
- Use a clean, soft cloth to clean the surface of the device and LCD. Do not use solvents or immerse the device into water or other liquids.
- Always keep the thermometer within the storage temperature range (- 20°C to 55°C or - 4°F to 131°F) and humidity range ($\leq 93\%$ non-condensing).
- It is recommended to store the thermometer in a dry location free from dust. Do not expose the thermometer to direct sunlight, high temperature/ humidity or any extreme environment, otherwise the function will be reduced.
- When the ambient temperature of the thermometer changes too much, such as moving the thermometer from one place of lower temperature to another place of higher temperature, allow the thermometer to remain in a room for 30 minutes where the temperature is between 15°C to 40°C.

DISPOSAL

- Used batteries should not be disposed of in the household rubbish. Used Batteries should be deposited at a collection point.
- At the end of its life, the appliance should not be disposed of in household rubbish.
- Enquire about the options for environment-friendly and appropriate disposal. Take local regulations into account.

WARRANTY

Our company warrants Non-Contact Infrared Body Thermometer at the time of its original purchase and for the subsequent time period of 2 years.

The warranty does not cover the followings:

- The device series number label is torn off or cannot be recognized.
- Damage to the device resulting from misconnection with other devices.
- Damage to the device resulting from accidents.
- Changes performed by users without the prior written authorization of the company.
- Batteries and packaging are not covered under warranty

Note:

1. If you have any problems with this device, such as setting up, maintaining or using, please contact with SERVICE PERSONNEL of FKA Brands Ltd. Don't open or repair the device by yourself.
2. Please report to FKA Brands Ltd if any unexpected operation or events occur.
3. The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the user instructions, periodic re-adjustment is not required. If at any time you question the accuracy of temperature measurements, please contact us.
4. The patient is an intended operator. The patient can measure and change battery. Under normal circumstances and maintain the device and its accessories according to the user manual.

EMC DECLARATION

1. This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS; This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
2. * Caution: Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
3. *Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!

4. * Caution: this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

| Guidance and manufacturer's declaration–electromagnetic emissions | |
|---|-------------------|
| The Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Thermometer should assure that it is used in such an environment. | |
| Emissions test | Compliance |
| RF emissions CISPR 11 | Group 1 |
| RF emissions CISPR 11 | Class B |
| Harmonic emissions IEC 61000-3-2 | Not applicable |
| Voltage fluctuations Flicker emissions IEC 61000-3-3 | Not applicable |

Guidance and manufacturer's declaration–electromagnetic immunity

The Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Thermometer should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level |
|--|--|--|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 8 kV contact $\pm 2, \pm 4, \pm 8, \pm 15$ kV air | ± 8 kV contact $\pm 2, \pm 4, \pm 8, \pm 15$, kV air |
| Electrical fast transient/burst IEC 61000-4-4 | The input a.c. power ports: ± 2 KV The input d.c. power ports: ± 2 KV Signal input/output ports: ± 1 KV | Not applicable |
| Surge IEC 61000-4-5 | Input power ports: $+0.5, +1.0$ KV Signal input/output: $+2.0$ KV | Not applicable |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0.5 cycles for $> 95\%$ (sync angle (degrees): 0, 45, 90, 135, 180, 225, 270, 315) 1 cycles for $> 95\%$ UT (sync angle (degrees): 0) 25 (50Hz)/30 (60Hz) cycles for 30% U T (sync angle (degrees): 0) | Not applicable |
| Voltage interruption IEC 61000-4-11 | 250 (50Hz)/300 (60Hz) cycles for $> 95\%$ UT (sync angle (degrees): 0) | |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m | 30 A/m |

NOTE UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration—electromagnetic immunity

The Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Thermometer should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level |
|-------------------------------|--|---|
| Conducted RF IEC 61000-4-6 | 3Vrms 150kHz to 80 MHz | Not applicable |
| Radiated RF IEC 61000-4-3 | Professional healthcare environment: 3 V/m Home healthcare environment 10 V/m 80MHz - 2700 MHz | Professional healthcare environment: 3 V/m Home healthcare environment 10 V/m 80MHz - 2700 MHz ((Ⓢ)) |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as base stations from radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can not be predicted theoretic call with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Thermometer is used exceeds the applicable RF compliance level above, the Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Thermometer.
- b. Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration - RF wireless communication equipment immunity

| Test frequency (Mhz) | Band a) (MHz) | Service a) | Modulation b) | Modulation b) (W) | Distance (m) | IMMUNITY TEST LEVEL (V/m) |
|----------------------|---------------|--|--------------------------------------|-------------------|--------------|---------------------------|
| 385 | 380-390 | TETRA 400 | Pulse modulation b) 18Hz | 1.8 | 0.3 | 27 |
| 450 | 430-470 | GMRS 460, FRS 460 | FM c) \pm 5kHz deviation 1kHz sine | 2 | 0.3 | 28 |
| 710 | 704-787 | LTE Band 13,17 | Pulse modulation b) 217Hz | 0.2 | 0.3 | 9 |
| 745 | | | | | | |
| 780 | | | | | | |
| 810 | 800-960 | GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5 | Pulse modulation b) 18 Hz | 2 | 0.3 | 28 |
| 870 | | | | | | |
| 930 | | | | | | |
| 1720 | 1700-1990 | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4,25; UMTS | Pulse modulation b) 217 Hz | 2 | 0.3 | 28 |
| 1845 | | | | | | |
| 1970 | | | | | | |
| 2450 | 2400-2570 | Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7 | Pulse modulation 217 Hz | 2 | 0.3 | 28 |
| 5240 | 5100-5800 | WLAN 802.11 a/n | Pulse modulation 217 Hz | 0.2 | 0.3 | 9 |
| 5500 | | | | | | |
| 5785 | | | | | | |

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50% duty cycle square wave signal.

c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case

BODY TEMPERATURE

- Body temperature varies from person to person and fluctuates during the course of the day. For this reason, it is suggested to know one's normal, healthy forehead temperature to correctly determine the temperature.
- Body temperature runs approximately from 35.5°C to 37.8°C (95.9°F-100°F). To determine if one has a fever, compare the temperature detected with the person's normal temperature. A rise over the reference body temperature of 1°C (1°F) or more is generally indication of fever.
- Different measurement sites (rectal, axillary, oral, frontal, auricular) will give different readings. Therefore it is wrong to compare the measurement taken from different sites.
- Here below are typical temperatures for adults, based on different measurement sites:

Rectal: 36.6°C to 38°C /97.9°F-100.4°F

Axillary: 34.7°C to 37.3°C /94.5°F-99.1°F

ASTM laboratory accuracy requirements in the display range of 37°C (98°F to 102°F) for IR thermometers is ± 0.2 ($\pm 0.4^\circ\text{F}$), whereas for mercury in-glass and electronic thermometers, the requirement per ASTM Standards E667-86 and E1112-86 is ± 0.1 ($\pm 0.2^\circ\text{F}$).

Caution: This infrared thermometer meets requirements established in ASTM Standard (E1965-98) Except of clause 5.2.2. It displays subject's temperature over a range of 34.0°C ~43.0°C. Full responsibility for the conformance of this product to the standard is assumed by (Hetaida Technology Co., Ltd. Add: Room 801 ,802, 803, 804, 901, 2# Building Scientific Research Center, Songhu Intelligent Valley, No.6 Minfu Road, Liaobu Town, Dongguan City, Guangdong Province, P.R.China).

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FKA Brands Ltd, Somerhill Business Park, Five Oak Green Road,
Tonbridge, Kent TN11 0GP.
Customer Service: service@homedics.co.uk

Share Info Consultant Service LLC Repräsentanbüro
Heerdter Lohweg 83, 40549 Düsseldorf

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2 YEAR GUARANTEE

FKA Brands Ltd guarantees this product from defect in material and workmanship for a period of 2 years from the date of purchase, except as noted below. This FKA Brands Ltd product guarantee does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorised accessory; alteration to the product; or any other conditions whatsoever that are beyond the control of FKA Brands Ltd. This guarantee is effective only if the product is purchased and operated in the UK / EU. A product that requires modification or adaptation to enable it to operate in any country other than the country for which it was designed, manufactured, approved and / or authorised, or repair of products damaged by these modifications is not covered under this guarantee. FKA Brands Ltd shall not be responsible for any type of incidental, consequential or special damages. To obtain guarantee service on your product, return the product post-paid to your local service centre along with your dated sales receipt (as proof of purchase). Upon receipt, FKA Brands Ltd will repair or replace, as appropriate, your product and return it to you, post-paid. Guarantee is solely through HoMedics Service Centre. Service of this product by anyone other than the HoMedics Service Centre voids the guarantee. This guarantee does not affect your statutory rights. For your local HoMedics Service Centre, go to www.homedics.co.uk/servicecentres

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Hetaida Technology Co., Ltd. 4F,
FKA Brands Ltd, 29 Earlsfort Terrace, Dublin 2, Ireland.